



intended to be used for moistening and decontamination by physical cleansing of the nasal vestibules as well as for supportive treatment of irritated skin underneath the nasal opening.

octenisan® md nasal gel

Our Plus

- decontamination by physical cleansing of the nasal vestibules
- to support treatment of irritated skin underneath the nasal opening
- for moistening
- excellent compatibility

CE 0499

Application areas

The hydrophilic properties of octenisan® md nasal gel ensure efficient moisturising of the nasal vestibules. The moisturising effects of hydrogel promote a favourable physiological milieu which prevents the settlement of bacteria. Furthermore, restoring the physiological functions of protection and cleaning contribute to effective decontamination of e.g. particle, skin residuals and multiresistant pathogens. Irritated skin underneath the nasal opening is also supplied with moisture through the hydrogel and thus can heal optimally. Epithelial cells must be able to proliferate and migrate for the healing process. This is best possible in a moist wound environment.

Instructions for use

Unless otherwise prescribed, use octenisan® md nasal gel up to twice a day as needed. octenisan® md nasal gel is for use in the nasal vestibules only. Apply the gel e.g. with a cotton swab in a sufficient quantity to the frontal area of the nose (nasal vestibules), so that the surface is completely covered. For this purpose insert octenisan® md nasal gel into the nostril, close the nose by pressing on the side and then massage between thumb and finger. This ensures an even distribution of the gel. If needed, the nose should be thoroughly cleaned before with a wet cotton swab. Excessive gel is removed with a cotton swab. Make sure that a moist film remains. Direct contact with the opening of the tube must be avoided to prevent contamination with germs.



octenisan[®] md nasal gel

Product data

Ingredients:

Purified water, Propylene Glycol, Hydroxyethylcellulose, Octenidine HCl

Special advice

For all of the related safety information please refer to the Safety Data Sheet. To view the Safety Data Sheet please visit www.schuelke.com

octenisan[®] md nasal gel must not be used with known or suspected allergy to one of the ingredients (propylene glycol, hydroxyethyl cellulose, octenidine HCl). Consult with your physician, if you have any doubts. Do not use in combination with cleansing soaps, ointments, oils, enzymes or similar, as this may affect the preservation. Do not use in combination with PVP iodine, as this can lead to discoloration and impaired antiseptic effect of the PVP iodine. Do not use the octenisan[®] md nasal gel without medical supervision for longer than 2 weeks without interruption. octenisan[®] md nasal gel must not be used on cartilage, in the eyes or in the ear to avoid possible tissue damage. Do not insert too deep in the nose. Do not use for infusion or injection! Do not swallow. Do not use if tube is damaged. Store the tube at room temperature protected from direct sunlight.

Information for order

Item	Delivery form	Item no.
octenisan md nasal gel EN/PL/PT 6 ml TB	20/Carton	70002424
octenisan md Nasengel EN/PL 6 ml TB	20/Carton	70003046
octenisan md Nasengel EN/GR/AR 6 ml TB	20/Carton	70002735

Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

Expert opinion and information

Please visit our website for an overview of all available literature/reports on all schulke products online at <http://www.schuelke.com/>

For any queries you might have about the product or using it please call 0114 2543500 or E-mail mail.uk@schuelke.com



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke UK
Schülke & Mayr UK Ltd
Cygnets House, Meadowhall
Sheffield, S9 1AT
Phone +44-114-254 35 00
Fax +44-114-254 35 01
mail.uk@schuelke.com

schülke Headquarters
Schülke & Mayr GmbH
Robert-Koch-Str. 2
22851 Norderstedt, Germany
Phone +49 40 - 52100 - 0
Fax +49 40 - 52100 - 318
www.schuelke.com
info@schuelke.com