



# gigasept® FF (new)

gigasept® FF (new) is a high-level disinfectant based on succinic dialdehyde for manual reprocessing of surgical and sensitive medical instruments.

## Our plus:

- formaldehyde and GDA-free
- broad microbiological efficacy
- outstanding material compatibility
- 14-days standing time

## Application areas

Disinfection of thermolabile and thermostable instruments, e.g. flexible and rigid endoscopes, waterproof ultrasound probes, anaesthesia equipment, breathing masks, instruments etc.

## Instructions for use

gigasept® FF (new) is a concentrate and is diluted with cold water to the desired working concentration.

### Application example:

1 litre of a 6 % working solution corresponds to 940 ml of water and 60 ml of gigasept® FF (new).

Immerse endoscopes and instruments to be reprocessed in the gigasept® FF (new) solution. Ensure complete covering of the surfaces, also of hollow instruments, and allow them to soak in the solution. After disinfection, thoroughly rinse all instruments using water of at least drinking water quality, or preferably sterile distilled water or fully deionised water, to completely remove any residues of the solution. Please follow the reprocessing recommendations of the instrument manufacturers.

### Standing time:

A standing time of 14 days applies if the working solution is only used for the disinfection of previously cleaned instruments. National regulations may limit the timespan during which the solution can be used.

Apart from the manual treatment of endoscopes and instruments, gigasept® FF (new) is also suitable for semiautomated machines that operate by a cyclic process at room temperature.

## Microbiological efficacy

Efficacy	Concentration	Contact time
<b>bactericidal</b> EN 13727, EN 14561 in accordance with VAH - clean conditions	5 % (50 ml/l) 4 % (40 ml/l) 2 % (20 ml/l)	15 min. 30 min. 60 min.
<b>yeasticidal</b> EN 13624, EN 14562 in accordance with VAH - clean conditions	5 % (50 ml/l) 4 % (40 ml/l) 2 % (20 ml/l)	15 min. 30 min. 60 min.
<b>tuberculocidal</b> EN 14348, EN 14563 in accordance with VAH - clean conditions	5 % (50 ml/l) 4 % (40 ml/l) 2 % (20 ml/l)	15 min. 30 min. 60 min.
<b>mycobactericidal</b> EN 14348, EN 14563 in accordance with VAH - clean conditions	7 % (70 ml/l) 5 % (50 ml/l) 3 % (30 ml/l)	15 min. 30 min. 60 min.
<b>virucidal against enveloped viruses</b> (incl. HIV, HBV, HCV) EN 14476, EN 17111 - clean conditions	1 % (10 ml/l)	15 min.
<b>virucidal against enveloped viruses</b> (incl. HIV, HBV, HCV) in accordance with DVV (German Registered Association for Combating Viral Diseases)/RKI Guideline	1 % (10 ml/l)	15 min.
<b>fungicidal</b> EN 13624, EN 14562 - clean conditions	6 % (60 ml/l)	60 min.
<b>virucidal</b> EN 14476, EN 17111 - clean conditions	6 % (60 ml/l)	60 min.
<b>virucidal</b> in accordance with DVV (German Registered Association for Combating Viral Diseases)/RKI Guideline	8 % (80 ml/l)	60 min.
<b>sporicidal (C. difficile)</b> EN 17126 - clean conditions	6 % (60 ml/l)	8 h

## Listings

- IHO listed
- VAH certificate

CE 0297



# gigasept® FF (new)

## Product data

### Composition:

100 g solution contain the following active substances:  
93.9 g Reaction product of DMO-THF, ethanol and water.

### Chemical-physical data

Density	ca. 1.01 g/cm <sup>3</sup> / 20 °C
Color	green
Flash point	38.5 °C / Method: DIN 51755 Part 1
Form	liquid
pH	6.3 – 6.6 / 100 % / 20 °C
Viscosity, dynamic	No data available

## Special advice

### Always read the label and product information before use.

For professional use only. Any serious incidents related to the product are to be reported to the manufacturer and competent authority.

### Material compatibility:

Metals, rubber, glass, plastics, flexible endoscopes.

### Interactions:

As a general rule, do not mix gigasept® FF (new) with other products. If gigasept® FF (new) comes into contact with application solutions containing amines, there may be precipitation and irreversible discolorations, for example if cleaning agents containing amines are used during precleaning and are not adequately rinsed off.

### Storage:

Do not store for long periods above 25 °C. Keep container tightly closed. Protect from frost, heat and sunlight. Store in the original container. Do not use the product after expiry date.

## Information for order

Item	Delivery form	Item no.
gigasept® FF (new) -INT- 2 l-bottles	5 / Carton	on request
gigasept® FF (new) -INT- 5 l-canister	1 / Canister	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Accessories

Item	Item no.
Can key for 5 + 10 l	135810
Measuring beaker 500 ml	136101
Measuring beaker 50 ml	136102

## Related products

- gigazyme®
- gigazyme® X-tra

## Environmental information

Schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining out high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

Email: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

### schülke Headquarters

Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt  
Germany  
Phone +49 40 52100-0  
Fax +49 40 52100-318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)